



CS CERAMIC CO.,LTD

FAX:0086-0731-23455639 TEL:0086-0731-23455639

ADD:Chashan Town Liling City Hunan Province ,China Email: info@csceramic.com

Web:www.csceramic.com

Platinum Crucible

Item NO.	Platinum %	Capacity (ml)	Top D(mm)	Bottom D (mm)	Height (mm)	Thickness(mm)	Crucible Weight (g)	Lid weight(g)	Picture
CS-PT-1	≥99.95%	5	20	17	15	0.22	6	2	
CS-PT-2	≥99.95%	10	28	17	28	0.26	11	2.5	
CS-PT-3	≥99.95%	15	32	19	32	0.26	15	3	
CS-PT-4	≥99.95%	20	35	21	35	0.27	20	3.5	
CS-PT-5	≥99.95%	25	37	23	37	0.27	25	4	
CS-PT-6	≥99.95%	30	39	26	40	0.29	30	4.5	
CS-PT-7	≥99.95%	35	40	26	40	0.29	35	5	
CS-PT-8	≥99.95%	40	42	26	42	0.29	40	5	
CS-PT-9	≥99.95%	50	44	26	44	0.31	50	5.5	
CS-PT-10	≥99.95%	60	49	29	49	0.33	60	7	
CS-PT-11	≥99.95%	80	54	37	54	0.35	70	8	
CS-PT-12	≥99.95%	100	60	37	60	0.37	85	10	
CS-PT-13	≥99.95%	150	70	41	70	0.38	156	12	
CS-PT-14	≥99.95%	200	74	43	74	0.38	175	15	
CS-PT-15	≥99.95%	250	78	53	78	0.39	198	30	
CS-PT-16	≥99.95%	300	86	57	83	0.4	240	36	

Platinum Crucible Tongs -Stainless steel tong+Platinum coated

Type	Stainless steel tong specs	Stainless steel tong	Platinum coated weight (g)	Picture
CS-ST-1	6"	15	4	
CS-ST-2	8"	20	5	
CS-ST-3	10"	26	6	
CS-ST-4	12"	30	7	
CS-ST-5	14"	36	8	
CS-ST-6	16"	41	9	
CS-ST-7	18"	46	10	
CS-ST-8	20"	51	11	
CS-ST-9	22"	56	12	